LabChip GX Touch 培训手册

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9. 做 RNA 样品检测时，还需 DEPC treated water (nuclease free) and PCR cap strips used to cover the micro-plates

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HT DNA 1K Specifications

	HT DNA 1K
Sizing Range	25 bp - 1000 bp
Sizing Resolution ¹	± 15% from 25 - 100 bp
	± 10% from 100 - 150 bp, 700 - 1000 bp
	± 5% from 150 - 700 bp
Sizing Accuracy	± 10 %
Sizing Precision	5% CV
Linear Concentration Range	0.1 ng/μL - 50 ng/μL per fragment
Sensitivity	0.1 ng/μL
Maximum Total DNA	80 ng/μL, 50 ng/μL per fragment
Carry-Over	< 0.25%
Quantification Accuracy	± 30% or ± 1 ng/μL, whichever is greater
Quantification Precision	20% from 25 - 500 bp, 10% from 500 - 1000 bp
Maximum Salt Concentration	125 mM
Additives	BSA/ detergents should not exceed 0.05mg/mL/ 0.01% v/v)
Chip Lifetime ²	2000 samples per chip
Reagent Kit Lifetime	10 chip preps for 384 samples each, or 20 chip preps for 48 samples
Samples per Chip Prep	384 samples for 10 chip preps, 48 samples for 20 chip preps
Standard Assay: Specifications are defined for this Assay	HT DNA 1K Standard: For sizing of DNA fragments in 25 to 1000 base pair range. (Analysis time per sample - 68 seconds)
Extra Assays	HT DNA 1K High Resolution: For sizing of DNA fragments in 25 to 1000 base pair range. Greater resolution with longer analysis time per sample. (Analysis time per sample - 120 seconds)

HT DNA 5K Specifications

	HT DNA 5K
Sizing Range	100 bp - 5000 bp
Sizing Resolution ¹	± 10% from 150 - 500 bp
	± 15% from 100 - 150 bp, 500 - 1500 bp
	± 20% from 1500 - 5000 bp
Sizing Accuracy	± 10 %
Sizing Precision	5% CV
Linear Concentration Range	0.25 ng/μL - 50 ng/μL per fragment
Sensitivity	0.25 ng/μL
Maximum Total DNA	80 ng/μL total, 50 ng/μL per fragment
Carry-Over	< 0.5%
Quantification Accuracy	± 30% or ± 1 ng/μL, whichever is greater
Quantification Precision	20% CV
Maximum Salt Concentration	125 mM
Additives	BSA/ detergents should not exceed 0.05mg/mL/ 0.01% v/v)
Chip Lifetime ²	2000 samples per chip
Reagent Kit Lifetime	10 chip preps for 384 samples each, or 20 chip preps for 48 samples
Samples per Chip Prep	384 samples for 10 chip preps, 48 samples for 20 chip preps
Standard Assay: Specifications are defined for this Assay	HT DNA 5K Standard: For sizing of DNA fragments in 100 to 5000 base pair range. Fastest analysis time per sample compared with all available assays. (Analysis time per sample - 28 seconds)

¹Resolution is defined as half height or better separation of two peaks. Actual separation performance can depend on the sample and application. Peaks that are resolved less than half height can still be accurately identified by the system software.

²Expected chip lifetime is based on use under normal laboratory conditions and adherence to Caliper preparation protocols, sample guidelines and storage conditions. Individual results may vary.

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HT DNA Reagent Kit

- DNA 1K, Part Number CLS760673
- DNA 5K, Part Number CLS760675

The kit contains enough reagents for 20 Low-throughput (LT) or 10 High-throughput (HT) chip preparations. Up to 48 samples can be tested with a LT chip preparation. Up to 384 samples can be tested with a HT chip preparation.

Reagents	Label	Quantity
----------	-------	----------

DNA Dye Concentrate	Blue	1 vial, 0.09 mL
Chip Storage Buffer	White	9 vials, 1.8 mL each
DNA Gel Matrix	Red	5 vials, 1.1 mL each
10X DNA Ladder	Yellow	1 vial, 0.26 mL
DNA Marker	Green	1 vial, 1.2 mL

Consumable Items	Supplier Number	Quantity
Spin Filters	Costar, Cat. # 8160	10
Ladder Tubes	Genemate, Cat. # C-3258-1	20, 0.2 mL
Detection Window Cleaning Cloth	VWR, Cat. # 21912-046	1
Swab	ITW Texwipe®, Cat. # TX758B	3
Buffer Tubes	E&K Scientific, Cat. # 697075-NC	20, 0.75 mL

HT DNA LabChip, Version 2

- DNA 1K Chip, Part Number 760517
- DNA 5K Chip, Part Number 760435

Item	Quantity
DNA Chip	1

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